## Transfer protocol of human HSC into NOG mice

#### 8-12 week old NOG mice irradiated with 2.5 Gy Intravenous transfer of 1 - 0.5 x 10<sup>5</sup> hCD34<sup>+</sup> cells Newborn NOG mice irradiated with 1 Gy I - 0.5 x 10<sup>5</sup> hCD34<sup>+</sup> cells Newborn NOG mice irradiated with 1 Gy I - 0.5 x 10<sup>5</sup> hCD34<sup>+</sup> cells I - 0.5 x

Intrahepatic or intravenous transfer of 1 - 0.1 x 10<sup>5</sup> hCD34<sup>+</sup> cells

#### Mice:

#### **SINCE 1952**

CIEA

Adult NOG mice are aged 8 - 12 weeks. Newborn mice are 1 – 2 days old.

#### Irradiation:

Irradiation is performed with 2 - 2.5 Gy to adult mice and with 1 Gy to newborns one day prior to cell transfer under SPF conditions. Mice less than 10 g body weight may die with this dose irradiation.

#### Cell preparation

- 1. When CD34+ cells (e.g. 1 x  $10^{6}$ /tube) are frozen, cells are thawed in a  $37^{\circ}$ C water bath.
- 2. Transfer cells into 50 mL conical tube and add 1 mL of PBS containing 2% fetal bovine serum (2%FBS-PBS) drop by drop shaking the tube slowly.
- 3. After further addition of 18 mL of 2%FBS-PBS, centrifuge it at 1,200 rpm for 5 min at room temperature.
- 4. Re-suspend the cells in 10 mL of 2%FBS-PBS, re-centrifuge.
- 5. Re-suspend in 2 mL of PBS.
- 6. Using  $15 \,\mu$ L of the suspension, count the cell number and viability after adding  $15 \,\mu$ L of 0.25% trypan-blue solution. Add a further 3 mL of PBS to reach a total of 5 mL.
- 7. Transfer 1 mL of the cell suspension into each 1.5 mL cryotube.
- 8. Put the tubes into a 50 mL conical tube and carry it to the animal facility.
- 9. After dipping it in an antiseptic solution, carry it into the animal room.

#### Cell transfer

#### For adult mice

1. 0.25 mL (1 - 0.5 x 10<sup>4</sup>) of the cell suspension is injected into mice via the tail vein with a 1 mL syringe with 27 G or a Microinjector syringe with 29 G under slight anesthesia with Isoflurane.

#### For newborn mice

- 1. Slightly anesthetize the newborn with Isoflurane.
- 2. Hold the newborn with left hand so the head of the mouse face right.
- 3.  $0.1 \text{ mL} (1-5 \text{ x } 10^4)$  of the cell suspension is injected via the face vein with a Microinjector syringe with 29 G.
- 4. Return the newborn to the mother gently, do not give it a shock.

## Engraftment of human cell in NOG mice



**SINCE 1952** 



8 - 20 weeks after intravenous transplantation of 5 x 10<sup>4</sup> human cord blood derived CD34+ cells into NOG and NOD-*scid* mice, human CD45+ hematopoietic cells in peripheral blood of mice were detected in flow cytometry. The high engraftment rate of human cells was observed in NOG mice.

Table 1. Limiting dilution assay in NOG mice.		
CD34+ cell dose	No. mice transplanted	No. mice successfully engrafted*
1,000	5	5
200	3	2
100	6	3

\* Successful engraftment was defined the presence of at least 0.1% human CD45+ cells in bone marrow by flow cytometory.

In order to evaluate the frequency of CB CD34<sup>+</sup> cells capable of engrafting NOG mice, these mice were transplanted with CB CD34<sup>+</sup> cells in a limiting dose. All recipient mice transplanted with 1 x 10<sup>3</sup> CB CD34<sup>+</sup> cells showed successful engraftment (>0.1%). Further, as few as 100 cells could be engrafted in 2 of 6 mice and their multi-lineage differentiation.

## Multi-lineage cell differentiation from HSCs in NOG mice

#### In bone marrow

CIEA SINCE 1952



The rate of human CD45<sup>+</sup> cells in bone marrow of transplanted mice and the high efficacy of multi-lineage cell differentiation. 11 weeks after CD34<sup>+</sup> cell transplantation, bone marrow were removed from NOD/SCID/ $\gamma_c^{null}$  and NOD/Shi-*scid* mice treated with anti-asialo GM1 antibody two days before transplantation, and subjected to flow cytometry. NOD/SCID/ $\gamma_c^{null}$  showed significantly higher percentage of CD45<sup>+</sup> cells. Multi-lineage cells have been differentiated from CD34<sup>+</sup> cells with high efficacy (\* vs \*\*: P<0.01).

- Ito, M., H. Hiramatsu, K. Kobayashi, K. Suzue, M. Kawahata, K. Hioki, Y. Ueyama, Y. Koyanagi, K. Sugamura, K. Tsuji, T. Heike, and T. Nakahata. 2002. NOD/SCID/gamma(c)(null) mouse: an excellent recipient mouse model for engraftment of human cells. *Blood* 100:3175-3182.
- Yahata, T., K. Ando, Y. Nakamura, Y. Ueyama, K. Shimamura, N. Tamaoki, S. Kato, and T. Hotta. 2002. Functional human T lymphocyte development from cord blood CD34+ cells in nonobese diabetic/Shi-scid, IL-2 receptor gamma null mice. *J Immunol* 169:204-209.
- 3. Saito, Y., Y. Kametani, K. Hozumi, N. Mochida, K. Ando, M. Ito, T. Nomura, Y. Tokuda, H. Makuuchi, T. Tajima, and S. Habu. 2002. The in vivo development of human T cells from CD34(+) cells in the murine thymic environment. *Int Immunol* 14:1113-1124.
- 4. Hiramatsu, H., R. Nishikomori, T. Heike, M. Ito, K. Kobayashi, K. Katamura, and T. Nakahata. 2003. Complete reconstitution of human lymphocytes from cord blood CD34+ cells using the NOD/SCID/gammacnull mice model. *Blood* 102:873-880.
- 5. Watanabe, Y., T. Takahashi, A. Okajima, M. Shiokawa, N. Ishii, I. Katano, R. Ito, M. Ito, M. Minegishi, N. Minegishi, S. Tsuchiya, and K. Sugamura. 2009. The analysis of the functions of human B and T cells in humanized NOD/shiscid/{gamma}cnull (NOG) mice (hu-HSC NOG mice). *Int Immunol.*

from Ito, M. et al. NOD/SCID/gamma(c)(null) mouse: an excellent recipient mouse model for engraftment of human cells. *Blood* 100, 3175-3182. (2002).

## Development of lymphoid lineage cells in NOG mice transferred with human CD34+ cells





The above figure shows human cells in mouse peripheral blood 20 weeks after intravenous human HSC cell transfer. Human T and B cells were developed in NOG mice. IgM+IgD+ B cells and CD4+ and CD8+ T cells were also differentiated in NOG mice.



become dominant in the periphery including peripheral blood and spleen.

from Watanabe, S. et al. Hematopoietic stem cell-engrafted NOD/SCID/IL2Rgamma null mice develop human lymphoid systems and induce long-lasting HIV-1 infection with specific humoral immune responses. Blood 109, 212-218 (2007).

from Dr. Habu S. Tokai Univ.

## Development of myeloid lineage cells in NOG mice transferred with human CD34+ cells

## CIEA SINCE 1952



*from* Watanabe, S. *et al.* Hematopoietic stem cell-engrafted NOD/SCID/IL2Rgamma null mice develop human lymphoid systems and induce long-lasting HIV-1 infection with specific humoral immune responses. *Blood* 109, 212-218 (2007).

### Megakaryocyte

#### In bone marrow



hCD42b staining H; human megakaryocytes M; mouse megakaryocytes

### Mast cell



hTryptase staining

### Platelet

#### In peripheral blood



from Dr. Miyakawa Y, Keio Univ.

## **Development of lymphoid follicle in NOG mice reconstituted human CD34+ cells**





#### HLA staining

Human lymphoid follicle like structures (White arrow) were observed in spleen of NOG mice 18 wks after human HSC transfer.

FDC-M1 murine FDC CD3 T cells **CD205 DC** 

*from* Watanabe, S. *et al.* Hematopoietic stem cell-engrafted NOD/SCID/IL2Rgamma null mice develop human lymphoid systems and induce long-lasting HIV-1 infection with specific humoral immune responses. *Blood* 109, 212-218 (2007).

Human B, T cells and other monocytes/macrophages and DCs, but not FDCs were associated with the generation of lymphoid follicle-like structures observed in spleen.

## Analysis of T cell receptor V $\beta$ repertoire



**SINCE 1952** 

 $V\beta$ Four to six months after transplantation of 2×10<sup>4</sup> to 5×10<sup>4</sup> CB CD34<sup>+</sup> cells, spleen cells were taken and the TCR Vß repertoire was analyzed by flow cytometry using a panel of 24 different antibodies. The results of three independent experiments are shown.

*from* Hiramatsu, H., R. Nishikomori, T. Heike, M. Ito, K. Kobayashi, K. Katamura, and T. Nakahata. 2003. Complete reconstitution of human lymphocytes from cord blood CD34+ cells using the NOD/SCID/gammacnull mice model. *Blood* 102:873-880.

## Activation of human T cells developed in NOG mice

CIEA SINCE 1952

## Proliferation and cytokine responses



## Activation markers



from Dr. Habu S, Tokai Univ.

## Antibody production in NOG mice transferred human HSCs

# CIEA



#### Non-specific antibodies

Anti-DNP-KLH specific antibodies





# Summary of human hematopoietic cells differentiated from CB CD34<sup>+</sup> cells in NOG mice



